


Pool grips – designed with midwives

-taking into account the hygienic regulations –



Specially developed open chromed handles are installed outside the water filling of the birthing and relaxation bathtub and comply with the highest hygienic regulations.

The open, fully encapsulated fixings of the handles are located outside the interior of the bathtub, thus avoiding permanent contamination with sewage.


An ergonomic positioning of the open chromed handles, gives the hands of the mother the possibility to form a line with the forearm. The body-friendly grip position allows to take an active birth posture.

Compared to the closed handles a simpler and more thorough cleaning of the open chrome handles is ensured.



The swivelling cord Variotrac® in combination with the birthing and relaxation pool

Equipment



Inflow tap
With 2 pluggable inflow arcs for a faster filling of the tub




PALL- Filter
For prevention of water germ infection and legionella infection.




Option for stair landing
2 castors at the handrail, by tilting the staircase manoeuvring is very easy



Under water spotlights
Provide subtle lighting for low-light areas



Lighting set spectrallight
Ensures a pleasant bathing experience through atmospheric lighting effects in the water



heating
Prevents a quick cooling of the bath water.



Integrated holding system
Rotatable, incl. 5 eyelets, carabiner and rope



PELVITUB®-cushion
Properly fixing by 9 suction cups, multifunctional in use



febromed

Febromed GmbH & Co. KG

Zum Mühlenteich 1

D-59302 Oelde

Phone +49 (0) 2522-30532

Fax +49 (0) 2522-308489

Web: www.febromed.com

E-Mail: info@febromed.de